

Product Datasheet

Human FOLH1 / PSMA, His Tag (orb1805368)

Description	Human FOLH1 / PSMA, His Tag is expressed from HEK293. It contains AA 707 (Accession No. Q04609).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	FOLH1 / PSMA
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Lys 44 - Ala 750
Isotype	His
Purity	> 85%
MW	80.62 kDa
Uniprot ID	Q04609
Expiration Date	6 months from date of receipt.