

## Product Datasheet

### Human FOLH1 / PSMA, Fc Tag (orb1805369)

<b>Description</b>	Human FOLH1 / PSMA, Fc Tag is expressed from HEK293. It contains AA 707 (Accession No. Q04609).
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/mg
<b>Tested Applications</b>	ELISA, Kinetics
<b>Target</b>	FOLH1 / PSMA
<b>Preservatives</b>	1×PBS, pH 6.0
<b>Storage</b>	-80°C
<b>Tag</b>	Fc
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Lys 44 - Ala 750
<b>Isotype</b>	Fc
<b>Purity</b>	> 85%
<b>MW</b>	211.2 kDa
<b>Uniprot ID</b>	<a href="https://www.uniprot.org/entry/Q04609">Q04609</a>
<b>Expiration Date</b>	6 months from date of receipt.