

Product Datasheet

Human PSCA, His Tag (orb1805370)

Description	Human PSCA, His Tag is expressed from HEK293. It contains AA 75 (Accession No. O43653).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	PSCA
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Leu 12 - Ser 86
Isotype	His
Purity	> 85%
MW	9.32 kDa
Uniprot ID	O43653
Expiration Date	6 months from date of receipt.