

Product Datasheet

Human MIF, His Tag (orb1805372)

Description	Human MIF, His Tag is expressed from HEK293. It contains AA 114 (Accession No. P14174).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	MIF
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Pro 2 - Ala 115
Isotype	His
Purity	> 85%
MW	7.43 kDa
Uniprot ID	P14174
Expiration Date	6 months from date of receipt.