

Product Datasheet

Human Integrin b7 / ITGB7, His Tag (orb1805374)

| | |
|----------------------------|---|
| Description | Human Integrin b7 / ITGB7, His Tag is expressed from HEK293. It contains AA 704 (Accession No. P2601... |
| Species/Host | Human |
| Reactivity | Human |
| Endotoxins | < 1 EU/mg |
| Tested Applications | ELISA, Kinetics |
| Target | Integrin b7 / ITGB7 |
| Preservatives | 1×PBS, pH 6.0 |
| Storage | -80°C |
| Tag | His |
| Note | For research use only |
| Protein Sequence | Glu 20 - His 723 |
| Isotype | His |
| Purity | > 85% |
| MW | 77.54 kDa |
| Uniprot ID | P26010 |
| Expiration Date | 6 months from date of receipt. |