

Product Datasheet

Human Integrin b7 / ITGB7, Fc Tag (orb1805375)

Description	Human Integrin b7 / ITGB7, Fc Tag is expressed from HEK293. It contains AA 704 (Accession No. P26010...
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	Integrin b7 / ITGB7
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	Fc
Note	For research use only
Protein Sequence	Glu 20 - His 723
Isotype	Fc
Purity	> 85%
MW	205.02 kDa
Uniprot ID	P26010
Expiration Date	6 months from date of receipt.