

Product Datasheet

Human GPC3 / Glypican-3, His Tag (orb1805376)

Description	Human GPC3 / Glypican-3, His Tag is expressed from HEK293. It contains AA 530 (Accession No. P51654)...
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	GPC3 / Glypican-3
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Gln 25 - Asn 554
Isotype	His
Purity	> 85%
MW	61.36 kDa
Uniprot ID	P51654
Expiration Date	6 months from date of receipt.