

Product Datasheet

Human TSPAN26 / CD37, His Tag (orb1805384)

Description	Human TSPAN26 / CD37, His Tag is expressed from HEK293. It contains AA 130 (Accession No. P11049).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	TSPAN26 / CD37
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Arg 112 - Asn 241
Isotype	His
Purity	> 85%
MW	15.99 kDa
Uniprot ID	P11049
Expiration Date	6 months from date of receipt.