

Product Datasheet

Human CD34, Fc Tag (orb1805387)

Description	Human CD34, Fc Tag is expressed from HEK293. It contains AA 259 (Accession No. P28906).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	CD34
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	Fc
Note	For research use only
Protein Sequence	Ser 32 - Thr 290
Isotype	Fc
Purity	> 85%
MW	107.12 kDa
Uniprot ID	P28906
Expiration Date	6 months from date of receipt.