

Product Datasheet

Human Neprilysin / CD10, His Tag (orb1805393)

Description	Human Neprilysin / CD10, His Tag is expressed from HEK293. It contains AA 699 (Accession No. P08473)...
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	Neprilysin / CD10
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Tyr 52 - Trp 750
Isotype	His
Purity	> 85%
MW	80.93 kDa
Uniprot ID	P08473
Expiration Date	6 months from date of receipt.