

## Product Datasheet

### Human ANGPTL3, His Tag (orb1805397)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human ANGPTL3, His Tag is expressed from HEK293. It contains AA 204 (Accession No. Q9Y5C1). |
| <b>Species/Host</b>        | Human   |
| <b>Reactivity</b>          | Human   |
| <b>Endotoxins</b>          | < 1 EU/mg   |
| <b>Tested Applications</b> | ELISA, Kinetics   |
| <b>Target</b>              | ANGPTL3   |
| <b>Preservatives</b>       | 1×PBS, pH 6.0   |
| <b>Storage</b>             | -80°C   |
| <b>Tag</b>                 | His   |
| <b>Note</b>                | For research use only   |
| <b>Protein Sequence</b>    | Ser 17 - Pro 220  |
| <b>Isotype</b>             | His   |
| <b>Purity</b>              | > 85%   |
| <b>MW</b>                  | 24.9 kDa  |
| <b>Uniprot ID</b>          | <a href="https://www.uniprot.org/entry/Q9Y5C1">Q9Y5C1</a>                                   |
| <b>Expiration Date</b>     | 6 months from date of receipt.  |