

Product Datasheet

Human Alpha-synuclein, His Tag (orb1805401)

Description	Human Alpha-synuclein, His Tag is expressed from HEK293. It contains AA 140 (Accession No. P37840).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	Alpha-synuclein
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Met 1 - Ala 140
Isotype	His
Purity	> 85%
MW	15.56 kDa
Uniprot ID	P37840
Expiration Date	6 months from date of receipt.