

Product Datasheet

Anti-TROP2 Reference Antibody (orb1805430)

Description Anti-TROP2 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the light C...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target TROP2

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only


Isotype IgG

Clonality Monoclonal


Purity > 95%

MW 146.06 kDa


Expiration Date 12 months from date of receipt.



Anti-TROP2
Reference
Antibody
(Sacituzumab...)



The purity of
Anti-TROP2
Reference
Antib...



Sacituzumab-
ADC
inhibited the
tumor growth...