



## Product Datasheet

### Anti-NCAM Reference Antibody (orb1805443)

<b>Description</b>	Anti-NCAM Reference Antibody and Monomethyl auristatin E (MMAE) used for the study of mul...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/mg
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FA, FACS, Kinetics
<b>Target</b>	NCAM
<b>Preservatives</b>	100 mM Pro-Ac 20 mM Arg pH 5.0
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95%
<b>MW</b>	145 kDa
<b>Expiration Date</b>	12 months from date of receipt.



Anti-NCAM1 / CD56 Reference Antibody (lo...



The purity of Anti-NCAM1 / CD56 Referenc...