



Product Datasheet

Anti-NCAM Reference Antibody (orb1805444)

Description	Anti-NCAM Reference Antibody and DM1, a potent derivative of maytansine, used for the stu...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	NCAM
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145 kDa
Expiration Date	12 months from date of receipt.



Anti-NCAM1 / CD56 Reference Antibody (lo...



The purity of Anti-NCAM1 / CD56 Referenc...