

---


## Product Datasheet

### Anti-ERBB2 / HER2 / CD340 Reference Antibody (orb1805445)

**Description**

Anti-ERBB2 / HER2 / CD340 Reference

the study of brea...



Anti-ERBB2 / HER2 / CD340 Reference Anti...

**Reactivity**

Human

**Endotoxins**

&lt; 1 EU/mg

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, FA, FACS, Kinetics

**Target**


ERBB2 / HER2 / CD340

**Preservatives**

100 mM Pro-Ac 20 mM Arg pH 5.0

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.



The purity of Anti-ERBB2 / HER2 / CD340 ...

**Note**

For research use only

**Isotype**

IgG

**Clonality**


Monoclonal

**Purity**

&gt; 95%

**MW**

145.48 kDa



Immobilized human HER2 His at 2 µg/mL c...

**Expiration Date**

12 months from date of receipt.