

## Product Datasheet

# Anti-ERBB2 / HER2 / CD340 Reference Antibody (orb1805446)

**Description** Anti-ERBB2 / HER2 / CD340 Reference AntibodyM1, a potent derivative of maytansine, used for the stud...

**Reactivity** Human

**Endotoxins** < 1 EU/mg

**Conjugation** Unconjugated

**Tested Applications** ELISA, FA, FACS, Kinetics

**Target** ERBB2 / HER2 / CD340

**Preservatives** 100 mM Pro-Ac 20 mM Arg pH 5.0

**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note** For research use only


**Isotype** IgG

**Clonality** Monoclonal


**Purity** > 95%

**MW** 145.48 kDa

**Expiration Date** 12 months from date of receipt.



Anti-ERBB2 / HER2 / CD340 Reference Anti...



The purity of Anti-ERBB2 / HER2 / CD340 ...