

Product Datasheet

Anti-TNFRSF8 / CD30 Reference Antibody (orb1805473)

Description

Anti-TNFRSF8 / CD30 Reference Antibody and


Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	TNFRSF8 / CD30
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.38 kDa
Expiration Date	12 months from date of receipt.



Anti-TNFRSF8 / CD30 Reference Antibody (...



The purity of Anti-TNFRSF8 / CD30 Refere...



Immobilized human CD30 His at 2 µg/mL c...