



Product Datasheet

Anti-LIV-1 / SLC39A6 Reference Antibody (orb1805475)


Description	Anti-LIV-1 / SLC39A6 Reference Antibody and Monomethyl auristatin E (MMAE) used for the study of brea...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	LIV-1 / SLC39A6
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-LIV-1 / SLC39A6 Reference Antibody ...



The purity of Anti-LIV-1 / SLC39A6 Refer...



Immobilized human LIV 1; His Tag at 2 µ...