
Product Datasheet

Anti-ERBB2 / HER2 / CD340 Reference Antibody (orb1805478)

Description

Anti-ERBB2 / HER2 / CD340 Reference

the study of H...

Reactivity

Human

Endotoxins

< 1 EU/mg

Conjugation

Unconjugated

Tested Applications

ELISA, FA, FACS, Kinetics

Target

ERBB2 / HER2 / CD340

Preservatives

100 mM Pro-Ac 20 mM Arg pH 5.0

Storage

-80°C

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

Purity


> 95%

MW


150 kDa

Expiration Date


12 months from date of receipt.



Anti-ERBB2 / HER2 / CD340 Reference Anti...



The purity of Anti-ERBB2 / HER2 / CD340 ...



Immobilized human HER2(23652) His at 2 ...