



Product Datasheet

Anti-SIRPa / CD172a Reference Antibody (orb1805487)

Description	Anti-SIRPa / CD172a Reference Antibody) is expressed from CHO. The heavy chain type is IgG4PE, and t...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	SIRPa / CD172a
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	144.36 kDa
Expiration Date	12 months from date of receipt.



Anti-SIRPa /
CD172a
Reference
Antibody
(...



The purity
of Anti-
SIRPa /
CD172a
Refere...