



Product Datasheet

Anti-TROP2 Reference Antibody (orb1805489)

Description	Anti-TROP2 Reference Antibody antibody and SN38 for the study of breast cancer.
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	TROP2
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-TROP2 Reference Antibody (Sacituzum...



The purity of Anti-TROP2 Reference Antib...