

Product Datasheet

Anti-TNFRSF10B / TRAILR2 / CD262 Reference Antibody (orb1805553)

Description Anti-TNFRSF10B / TRAILR2 / CD262 Reference Antibody is expressed from CHO. The heavy chain type is h...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target TNFRSF10B / TRAILR2 / CD262

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only


Isotype IgG

Clonality Monoclonal

Purity > 95%

MW 145.5 kDa

Expiration Date 12 months from date of receipt.



Anti-TNFRSF10B / TRAILR2 / CD262 Referen...



The purity of Anti-TNFRSF10B / TRAILR2 / ...