

## Product Datasheet

# Anti-TNFRSF10B / TRAILR2 / CD262 Reference Antibody (orb1805555)

**Description** Anti-TNFRSF10B / TRAILR2 / CD262 Reference Antibody is expressed from CHO. The heavy chain type is h...

**Reactivity** Human

**Endotoxins** < 1 EU/mg

**Conjugation** Unconjugated

**Tested Applications** ELISA, FA, FACS, Kinetics

**Target** TNFRSF10B / TRAILR2 / CD262

**Preservatives** 100 mM Pro-Ac 20 mM Arg pH 5.0

**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note** For research use only


**Isotype** IgG

**Clonality** Monoclonal

**Purity** > 95%

**MW** 145.5 kDa

**Expiration Date** 12 months from date of receipt.



Anti-TNFRSF10B / TRAILR2 / CD262 Referen...



The purity of Anti-TNFRSF10B / TRAILR2 / ...