
Product Datasheet

Anti-GDF15 / MIC1 Reference Antibody (orb1805630)

Description

Anti-GDF15 / MIC1 Reference Antibody is expressed

...

Reactivity

Human

Endotoxins

< 1 EU/mg

Conjugation

Unconjugated

Tested Applications

ELISA, FA, FACS, Kinetics

Target

GDF15 / MIC1

Preservatives

100 mM Pro-Ac 20 mM Arg pH 5.0

Storage

 Maintain refrigerated at 2-8°C for up to 2 weeks.
 For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

Purity


> 95%

MW


145.5 kDa

Expiration Date


12 months from date of receipt.



Anti-GDF15 / MIC1 Reference Antibody (po...



The purity of Anti-GDF15 / MIC1 Referenc...



Immobilized human GDF15 His at 2 µg/mL ...