

Product Datasheet

Anti-IL-23 Reference Antibody (orb1805698)

Description Anti-IL-23 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the light C...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target IL-23

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only


Isotype IgG

Clonality Monoclonal


Purity > 95%

MW 145.5 kDa

Expiration Date 12 months from date of receipt.



Anti-IL-23
Reference
Antibody
(LY2525623...



The purity of
Anti-IL-23
Reference
Antib...