



Product Datasheet

Anti-AXL / UFO Reference Antibody (orb1805842)


Description	Anti-AXL / UFO Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the lig...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	AXL / UFO
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.42 kDa
Expiration Date	12 months from date of receipt.



Anti-AXL / UFO Reference Antibody (enapo...



The purity of Anti-AXL / UFO Reference A...



Immobilized human Axl FC at 2 µg/mL can...