

Product Datasheet

Anti-CEACAM6 / CD66c Reference Antibody (orb1805848)


Description

Anti-CEACAM6 / CD66c Reference Antibody is


Reactivity	hulgG2, and t...
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	CEACAM6 / CD66c
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-CEACAM6 / CD66c Reference Antibody ...



The purity of Anti-CEACAM6 / CD66c Refer...



Immobilized human CEACAM6 His at 2 µg/m...