

---

## Product Datasheet

### Anti-TNFRSF4 / OX40 / CD134 Reference Antibody (orb1805884)

**Description**

Anti-TNFRSF4 / OX40 / CD134 Reference

chain type is hulg1...

**Reactivity**

Human

**Endotoxins**

&lt; 1 EU/mg

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, FA, FACS, Kinetics

**Target**

TNFRSF4 / OX40 / CD134

**Preservatives**

100 mM Pro-Ac 20 mM Arg pH 5.0

**Storage**

-80°C

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Monoclonal

**Purity**

&gt; 95%


**MW**

145.5 kDa


**Expiration Date**

12 months from date of receipt.


kDa NR M R



Anti-TNFRSF4 / OX40 / CD134 Reference An...



The purity of Anti-TNFRSF4 / OX40 / CD13...



Immobilized human TNFRSF4 His at 2 µg/m...