

Product Datasheet

Anti-B7-H4 / VTCN1 Reference Antibody (orb1805923)

Description Anti-B7-H4 / VTCN1 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target B7-H4 / VTCN1

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only


Isotype IgG

Clonality Monoclonal

Purity > 95%

MW 104.42 kDa


Expiration Date 12 months from date of receipt.



Anti-B7-H4 / VTCN1 Reference Antibody (a...



The purity of Anti-B7-H4 / VTCN1 Referen...



Immobilized human B7 H4 His at 2 µg/mL ...