



## Product Datasheet

# Anti-B7-H2 / ICOSL / CD275 Reference Antibody (orb1805925)


<b>Description</b>	Anti-B7-H2 / ICOSL / CD275 Reference Antibody is expressed from CHO. The heavy chain type is hulgG2,...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/mg
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FA, FACS, Kinetics
<b>Target</b>	B7-H2 / ICOSL / CD275
<b>Preservatives</b>	100 mM Pro-Ac 20 mM Arg pH 5.0
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95%
<b>MW</b>	145.5 kDa
<b>Expiration Date</b>	12 months from date of receipt.



Anti-B7-H2 / ICOSL / CD275 Reference Ant...



The purity of Anti-B7-H2 / ICOSL / CD275...



Immobilized human B7 H2 His at 2 µg/mL ...