


Product Datasheet

Anti-FcRn (FCGRT & B2M) Reference Antibody (orb1805934)


Description	Anti-FcRn (FCGRT & B2M) Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, an...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	FcRn (FCGRT & B2M)
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-FcRn (FCGRT & B2M) Reference Antibody



The purity of Anti-FcRn (FCGRT & B2M) Reference Antibody is > 95%.



Immobilized human b2M L FcRn at 2 µg/mL