

---


## Product Datasheet

### Anti-MSPR / RON / CD136 Reference Antibody (orb1805938)

**Description**

Anti-MSPR / RON / CD136 Reference Antibody is

kDa NR M R


 Anti-MSPR /  
RON /  
CD136  
Reference  
Antibo...

**Reactivity**

huIgG1, an...

Human

**Endotoxins**

&lt; 1 EU/mg

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, FA, FACS, Kinetics

**Target**

MSPR / RON / CD136

**Preservatives**

100 mM Pro-Ac 20 mM Arg pH 5.0

**Storage**

 Maintain refrigerated at 2-8°C for up to 2 weeks.  
For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Monoclonal

**Purity**



&gt; 95%

**MW**

145.5 kDa

**Expiration Date**

12 months from date of receipt.


 The purity  
of Anti-  
MSPR /  
RON /  
CD136  
Re...

 Hs746T  
cells were  
stained  
with Anti-  
MSPR...