

Product Datasheet

Anti-Complement C5aR1 Reference Antibody (orb1805939)


Description

Anti-Complement C5aR1 Reference Antibody is


kDa NR M R



Anti-Complement C5aR1 Reference Antibody...



The purity of Anti-Complement C5aR1 Refe...



Immobilized human C5AR1 VLP at 16 µg/mL

...

huIgG1, and ...

Reactivity

Human

Endotoxins

< 1 EU/mg

Conjugation

Unconjugated

Tested Applications

ELISA, FA, FACS, Kinetics

Target

Complement C5aR1

Preservatives

100 mM Pro-Ac 20 mM Arg pH 5.0

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

Purity

> 95%

MW

145.5 kDa

Expiration Date

12 months from date of receipt.