

---

## Product Datasheet


### Anti-TNFRSF10B / TRAILR2 / CD262 Reference Antibody (orb1805962)

**Description**


Anti-TNFRSF10B / TRAILR2 / CD262 Reference

|                            |                                 |
|----------------------------|---------------------------------|
| <b>Reactivity</b>          | Human                           |
| <b>Endotoxins</b>          | < 1 EU/mg                       |
| <b>Conjugation</b>         | Unconjugated                    |
| <b>Tested Applications</b> | ELISA, FA, FACS, Kinetics       |
| <b>Target</b>              | TNFRSF10B / TRAILR2 / CD262     |
| <b>Preservatives</b>       | 100 mM Pro-Ac 20 mM Arg pH 5.0  |
| <b>Storage</b>             | -80°C                           |
| <b>Note</b>                | For research use only           |
| <b>Isotype</b>             | IgG                             |
| <b>Clonality</b>           | Monoclonal                      |
| <b>Purity</b>              | > 95%                           |
| <b>MW</b>                  | 145.5 kDa                       |
| <b>Expiration Date</b>     | 12 months from date of receipt. |


kDa NR M R



Anti-TNFRSF10B / TRAILR2 / CD262 Referen...



The purity of Anti-TNFRSF10B / TRAILR2 /...



Immobilized Cynomolgus TRAIL R2/DR5/ TNF...