

Product Datasheet

Anti-TNFSF7 / CD27L / CD70 Reference Antibody (orb1805986)


Description

Anti-TNFSF7 / CD27L / CD70 Reference Antibody


Reactivity	hulgG1,...
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	TNFSF7 / CD27L / CD70
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-TNFSF7 / CD27L / CD70 Reference Ant...



The purity of Anti-TNFSF7 / CD27L / CD70...



Immobilized human CD70 FC at 4 µg/mL ca...