

---

## Product Datasheet




### Anti- FcRH5 / IRTA2 / CD307e Reference Antibody (orb1805993)

**Description**

Anti- FcRH5 / IRTA2 / CD307e Reference

<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/mg
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FA, FACS, Kinetics
<b>Target</b>	FcRH5 / IRTA2 / CD307e
<b>Preservatives</b>	100 mM Pro-Ac 20 mM Arg pH 5.0
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95%
<b>MW</b>	145.5 kDa
<b>Expiration Date</b>	12 months from date of receipt.

kDa NR M R


 Anti- FcRH5  
/ IRTA2 /  
CD307e  
Reference  
A...

 The purity  
of Anti-  
FcRH5 /  
IRTA2 /  
CD30...

 Immobilized  
human  
IRTA2 His  
at 2 µg/mL  
...