


Product Datasheet

Human TGF-Beta 2 ELISA Kit (orb180605)

Description	Human TGF-Beta 2 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A303-S414
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is 0-115 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	47748 MW
Uniprot ID	P61812
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human TGF-beta 2
PicoKine
ELISA Kit
stan...