

---

## Product Datasheet

### Anti-CSF2Rb / CD131 Reference Antibody (orb1806062)

**Description**

Anti-CSF2Rb / CD131 Reference Antibody is

huIgG4, and th...

**Reactivity**

Human

**Endotoxins**

&lt; 1 EU/mg

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, FA, FACS, Kinetics

**Target**

CSF2Rb / CD131

**Preservatives**

100 mM Pro-Ac 20 mM Arg pH 5.0

**Storage**

 Maintain refrigerated at 2-8°C for up to 2 weeks.  
 For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Monoclonal

**Purity**

&gt; 95%




**MW**

145.5 kDa

**Expiration Date**

12 months from date of receipt.

kDa NR M R


 Anti-  
 CSF2Rb /  
 CD131  
 Reference  
 Antibody  
 (...

 The purity  
 of Anti-  
 CSF2Rb /  
 CD131  
 Refere...

 Immobilized  
 human  
 CSF2RB His  
 at 2  
 µg/mL...