
Product Datasheet


Anti-CALCRL / CGRPR Reference Antibody (orb1806080)

Description


Anti-CALCRL / CGRPR Reference Antibody is

	hulgG2, and th...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	CALCRL / CGRPR
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.


kDa NR M R



Anti-CALCRL / CGRPR Reference Antibody (...



The purity of Anti-CALCRL / CGRPR Refere...



Immobilized human CGRPR VLP Protein at 1...