
Product Datasheet

Anti-SIRPa / CD172a Reference Antibody (orb1806109)

Description

Anti-SIRPa / CD172a Reference Antibody is

huIgG4SP, and ...

Reactivity

Human

Endotoxins

< 1 EU/mg

Conjugation

Unconjugated

Tested Applications

ELISA, FA, FACS, Kinetics

Target

SIRPa / CD172a

Preservatives

100 mM Pro-Ac 20 mM Arg pH 5.0

Storage

-80°C

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

Purity

> 95%




MW

145.5 kDa

Expiration Date

12 months from date of receipt.

kDa N M R


 Anti-SIRPa /
CD172a
Reference
Antibody
(...

 The purity
of Anti-
SIRPa /
CD172a
Refere...

 Immobilized
human
SIRPα V1
FC at4
µg/m...