
Product Datasheet




Anti-PDGFR α / CD140a Reference Antibody (orb1806125)

Description

 Anti-PDGFR α / CD140a Reference Antibody is

Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	PDGFR α / CD140a
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.

kDa N M R


 Anti-PDGFR α / CD140a Reference Antibody ...

 The purity of Anti-PDGFR α / CD140a Refer...

 Immobilized human PDGFR α / CD140a; Fc Ta...