



Product Datasheet

Anti-NKG2A / CD159a Reference Antibody (orb1806134)


Description	Anti-NKG2A / CD159a Reference Antibody is expressed from CHO. The heavy chain type is hulgG4SP, and ...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	NKG2A / CD159a
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-NKG2A /
CD94
Reference
Antibody
(mo...



The purity of
Anti-NKG2A /
CD94
Referenc...



Human
NKG2A/CD94
HEK293 cells
were stain...