
Product Datasheet

Anti-IL-7Ra / CD127 Reference Antibody (orb1806161)

Description

Anti-IL-7Ra / CD127 Reference Antibody is

huIgG4SP, and ...

Reactivity

Human

Endotoxins

< 1 EU/mg

Conjugation

Unconjugated

Tested Applications

ELISA, FA, FACS, Kinetics

Target

IL-7Ra / CD127

Preservatives

100 mM Pro-Ac 20 mM Arg pH 5.0

Storage

 Maintain refrigerated at 2-8°C for up to 2 weeks.
 For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

Purity


> 95%

MW

145.5 kDa

Expiration Date


12 months from date of receipt.



Anti-IL-7Ra / CD127 Reference Antibody (...



The purity of Anti-IL-7Ra / CD127 Refere...



Immobilized human IL-7Ra / CD127 His at ...