
Product Datasheet

Anti-IL-4Ra / CD124 Reference Antibody (orb1806164)

Description

Anti-IL-4Ra / CD124 Reference Antibody is

huIgG4SP, and ...

Reactivity

Human

Endotoxins

< 1 EU/mg

Conjugation

Unconjugated

Tested Applications

ELISA, FA, FACS, Kinetics

Target

IL-4Ra / CD124

Preservatives

100 mM Pro-Ac 20 mM Arg pH 5.0

Storage

 Maintain refrigerated at 2-8°C for up to 2 weeks.
 For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

Purity

> 95%


MW

145.5 kDa


Expiration Date

12 months from date of receipt.


kDa NR M R



Anti-IL-4Ra / CD124 Reference Antibody (...



The purity of Anti-IL-4Ra / CD124 Refere...



Immobilized human IL4RA His at 2 µg/mL ...