



Product Datasheet

Anti-IL-31Ra Reference Antibody (orb1806167)


Description	Anti-IL-31Ra Reference Antibody is expressed from CHO. The heavy chain type is hulgG2, and the light...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	IL-31Ra
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-IL-31Ra Reference Antibody (nemoliz...



The purity of Anti-IL-31Ra Reference Ant...



Immobilized human IL31R His at 2 µg/mL ...