


Product Datasheet

Mouse Relaxin 1 ELISA Kit (orb180617)

Description	Mouse Relaxin 1 ELISA Kit
Reactivity	Mouse
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: R23-R54 (B) & E161-C185 (A)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 78 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
MW	20571 MW
Uniprot ID	P47932
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse Relaxin 1 PicoKine ELISA Kit stand...