

Product Datasheet

Anti-IGF1R / CD221 Reference Antibody (orb1806181)

Description

Anti-IGF1R / CD221 Reference Antibody is

hulgG1, and the...



Anti-IGF1R / CD221 Reference Antibody (t...

Reactivity

Human

Endotoxins

< 1 EU/mg

Conjugation

Unconjugated

Tested Applications

ELISA, FA, FACS, Kinetics

Target


IGF1R / CD221

Preservatives

100 mM Pro-Ac 20 mM Arg pH 5.0

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.



The purity of Anti-IGF1R / CD221 Referen...

Note

For research use only

Isotype

IgG

Clonality


Monoclonal

Purity

> 95%

MW

145.5 kDa



Immobilized human IGF1R His at 2 µg/mL ...

Expiration Date

12 months from date of receipt.