


Product Datasheet

Human TFF2/Trefoil factor 2 ELISA Kit (orb180619)

Description	Human TFF2/Trefoil factor 2 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E24-Y129
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:4, concentration in serum and plasma is around 2000-3800 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	14284 MW
Uniprot ID	Q03403
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
TFF2
PicoKine
ELISA Kit
standard
C...