



Product Datasheet

Anti-GCGR Reference Antibody (orb1806207)


Description	Anti-GCGR Reference Antibody is expressed from CHO. The heavy chain type is hulg4SP, and the light ...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	GCGR
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.



Anti-GCGR Reference Antibody (crotedumab...



The purity of Anti-GCGR Reference Antibo...



Human GCCR His EGFP HEK293 cells were st...